


ELEVATORS & LUXATORS




E



ELEVATORS

Elevators are used in a variety of dental procedures, including routine extractions, removing impacted teeth, and loosening teeth that are malpositioned or at risk of fracture. They work on the principle of leverage to dislodge a tooth from its socket.

(SHOWN ACTUAL SIZE)



3701 HEIDBRINK 1



3702 HEIDBRINK 2 3704 SELDIN 12M 3705 SELDIN 34 3706 SELDIN 34S 3706S SELDIN 34S SERRATED 3707 HEIDBRINK 40 3708 HEIDBRINK 41 3709 APEXO 80



3710 APEXO 81 3711 APEXO 81E 3712 APEXO 301 3712S APEXO 301 SERRATED 3712W APEXO 301 WIDE 3712WS APEXO 301 WIDE SERRATED 3713 APEXO 302 CVD INNER 2.8mm ±5% 3714 SELDIN 2S



3715 APEXO 303 CVD INNER 2.8mm ±5% 3716 CRANE PICK 8 3717 SELDIN 1L 3718 SELDIN 1R 3719 SELDIN 4L 3720 SELDIN 4R 3721 CRYER M25 3722 CRYER M26



3725 CRYER S30A 3726 CRYER S31A 3729 CRYER L46L 3730 CRYER L46R 3731 BARRY 50 3732 BARRY 51 3734 MILLER APEXO 71 3735 MILLER APEXO 72

SHOWN ACTUAL SIZE



3745 FLOHR R-SHARP

SHOWN SMALLER THAN ACTUAL SIZE



3747S BERNARD 3 LEFT SMALL HANDLE

3750S BERNARD 4 RIGHT SMALL HANDLE



3736
MILLER
APEXO 73



3737
MILLER
APEXO 70



3738
MILLER
APEXO 74



3739
APICAL 2S/301
EXTRA THIN



3740
APICAL 2S/302
EXTRA THIN



3741
APICAL 4R/303
EXTRA THIN



3742
FLOHR
L-SHARP



3743
FLOHR
L-SHARP



3747
BERNARD 3
LEFT



3748
BERNARD 2
STR



3749
BERNARD 2
CVD



3749IN
BERNARD
IN



3749OUT
BERNARD
OUT



3750
BERNARD 4
RIGHT



3750D
FLAP
OPENER



3750DT
FLAP
OPENER
TITANIUM GOLD



3759
COPELAND 3
(4.5mm)



3760
COPELAND 2
(3.5mm)



3765
LINDO
STR 3mm



3766
LINDO
STR 4mm



3767S
SELDIN
46S IN
SERRATED



3768S
SELDIN
46S OUT
SERRATED



3787
FRIEDMAN
GOUGE
F31 OUT
SERRATED



3788
FRIEDMAN
GOUGE
F31 IN
SERRATED



3801
FLOHR
3R



3802
FLOHR
1S



3803
FLOHR
2L



3805
KOPP
1S



3806
KOPP
3R



3807
WWJL
SMALL HANDLE



3808
WWJS
SMALL HANDLE



3809
WWJR
SMALL HANDLE



3814
CRYER
39



3814S
CRYER 39
SMALL HANDLE



3815
CRYER
40R



3815S
CRYER 40R
SMALL HANDLE



3818S
E77R OUT
4mm
SERRATION



3820S
E77 INNER
4mm
SERRATION



3822
HYLIN
2S



3823
HYLIN
1S



3824
HYLIN
3S



3827 COLEMAN 2LS



3828
COLEMAN
2RS



3830
LINDO LEVIN
INNER
SERRATED



3831
LINDO LEVIN
OUTER
SERRATED



3841
COGSWELL
B



3842
WING 1
EXTRA THIN



3843
WING 2
EXTRA THIN



3844
WING
3



3845
WING
4



3852
POTTS
6



3853
POTTS
7



3837 BARRY 17L T-BAR



3838 BARRY 17R T-BAR



3839 WINTER 14L T-BAR



3840 WINTER 14R T-BAR



3850 HOURIGAN 4 T-BAR



3851 HOURIGAN 3 T-BAR



3852T POTTS 6 T-BAR



3853T POTTS 7 T-BAR



3856
ASL
2.5mm L



3857
ASL
2.5mm R



3878
HYLIN
12cm
EXTRA THIN



3878L
HYLIN
16cm
EXTRA THIN




3882 LIZA L




3883 LIZA R


4034
POWER FLEX ELEVATOR SET




4034-1 POWER FLEX APEXO 301 STR 2.8mm ±5% TITANIUM



4034-2 POWER FLEX APEXO 302 CVD IN 2.8mm ±5% TITANIUM




4034-3 POWER FLEX APEXO 303 CVD OUT 2.8mm ±5% TITANIUM



4034-4 POWER FLEX SELDIN 34S 4.0mm ±5% TITANIUM



4034-5 POWER FLEX E77 UP 4.0mm ±5% TITANIUM



4034-6 POWER FLEX E77 DOWN 4.0mm ±5% TITANIUM

4031


POWER TWIST ELEVATOR SET




4023 POWER TWIST ELEVATOR STR-L



4024 POWER TWIST ELEVATOR STR-R



4026 POWER TWIST ELEVATOR CVD-L



4027 POWER TWIST ELEVATOR CVD-R

3799

POWER FLEX PDL TITANIUM ELEVATORS SET OF 7



3799-1 POWER FLEX PDL TITANIUM ELEVATOR 3mm STR



3799-2 POWER FLEX PDL TITANIUM ELEVATOR 3mm CVD



3799-3 POWER FLEX PDL TITANIUM ELEVATOR SPADE



3799-4 POWER FLEX PDL TITANIUM SMALL DISTAL



3799-5 POWER FLEX PDL TITANIUM SMALL MESIAL




3799-6 POWER FLEX PDL TITANIUM LARGE DISTAL



3799-7 POWER FLEX PDL TITANIUM LARGE MESIAL

LUXATORS

Luxators are used to loosen and remove teeth from the mouth. They are sharp with a thin blade that's inserted into the periodontal ligament to cut the fibers and separate the tooth from the socket.



3751 2mm STR



3752
3mm STR



3753
5mm STR



3754
2mm CVD



3755
3mm CVD



3756
5mm CVD



3776
LUX-TOME
3mm



3777
LUX-TOME
4mm



3778
LUX-TOME
5mm



3779
3mm STR
SMALL HANDLE



3780
4mm STR
SMALL HANDLE



3781
5mm STR
SMALL HANDLE



3782
3mm CVD
SMALL HANDLE



3783
4mm CVD
SMALL HANDLE



3784
5mm CVD
SMALL HANDLE




3785
4mm STR
SMALL HANDLE




3786
4mm CVD
SMALL HANDLE


4035 POWER FLEX LUXATOR SET




4035-1 POWER FLEX LUXATOR STR 2.5mm ±5% TITANIUM



4035-2 POWER FLEX LUXATOR CVD 2.5mm ±5% TITANIUM




4035-3 POWER FLEX LUXATOR STR 3.5mm ±5% TITANIUM



4035-4 POWER FLEX LUXATOR CVD 3.5mm ±5% TITANIUM



4035-5 POWER FLEX LUXATOR STR 4.5mm ±5% TITANIUM



4035-6 POWER FLEX LUXATOR CVD 4.5mm ±5% TITANIUM

LUX-TOMES

When needing the ability to create space and luxate, but also have enough thinness to insert deep into the periodontal pockets, our Power Lux-Tomes offer a hybrid approach to performing extractions.



3769 ACCU LUX-TOME STR SERRATED TITANIUM



3770 ACCU LUX-TOME CVD SERRATED TITANIUM



3771 ACCU LUX-TOME SPADE SERRATED TITANIUM



3772 ACCU LUX-TOME DISTAL SERRATED TITANIUM



3773 ACCU LUX-TOME MESIAL SERRATED TITANIUM



3774 ACCU LUX-TOME LARGE DISTAL SERRATED TITANIUM



3775 ACCU LUX-TOME LARGE MESIAL SERRATED TITANIUM




3790 ACCU BERNARD SPADE OFFSET IN



3791 ACCU BERNARD SPADE OFFSET OUT



3792 ACCU BERNARD SPADE LEFT



3793 ACCU BERNARD SPADE RIGHT



3794 ACCU LUX-TOME T-BAR

POWER FLEX LUX-TOMES

The Power Lux-Tomes that you know and love, except enhanced with NiTi tips for flexibility, durability, and a lightweight handle design that allows for maximum control and usage inside the mouth.

3796 POWER FLEX SPADE KIT



3796-1 POWER FLEX SPADE STR



3796-2 POWER FLEX SPADE CVD



3796-3 POWER FLEX SPADE LEFT



3796-4 POWER FLEX SPADE RIGHT



3797-1 POWER FLEX ACCU LUX-TOME 3mm STR SERRATED



3797-2 POWER FLEX ACCU LUX-TOME 3mm CVD SERRATED



3797-4 POWER FLEX ACCU LUX-TOME SMALL DISTAL SERRATED



3797-5 POWER FLEX ACCU LUX-TOME SMALL MESIAL SERRATED



3797-6 POWER FLEX ACCU LUX-TOME MEDIUM DISTAL SERRATED



3797-7 POWER FLEX ACCU LUX-TOME MEDIUM MESIAL SERRATED

3798 POWER FLEX PDL ACCU LUXATING FULL SET OF 10

3796-1 POWER FLEX SPADE STR

3796-2 POWER FLEX SPADE CVD

3796-3 POWER FLEX SPADE LEFT

3796-4 POWER FLEX SPADE RIGHT

3797-1 POWER FLEX ACCU LUX-TOME 3mm STR SERRATED

3797-2 POWER FLEX ACCU LUX-TOME 3mm CVD SERRATED

3797-4 POWER FLEX ACCU LUX-TOME SMALL DISTAL SERRATED

3797-5 POWER FLEX ACCU LUX-TOME SMALL MESIAL SERRATED

3797-6 POWER FLEX ACCU LUX-TOME MEDIUM DISTAL SERRATED

3797-7 POWER FLEX ACCU LUX-TOME MEDIUM MESIAL SERRATED

ROOT TIP PICKS

Root tip picks are used to retrieve root tips that are in difficult to reach areas or that have broken away from the tooth during extraction. The sharp tip is placed into the periodontal ligament space and the root is teased from the socket.




3901 HEIDBRINK 1



3902 HEIDBRINK 2 3903 HEIDBRINK 3 3904 HEIDBRINK 4 3905 HEIDBRINK 5 3907 APICAL 9 3908 APICAL 9R 3916 HOWARD 1



3917 HOWARD 2 3918 HOWARD 3 3919 HOWARD 4 3920 HOWARD 5 4020 78 STR 4021 79L 4022 80R



3908N APICAL 9R 6EZ



3909E APICAL 9L 6EZ



3906 HEIDBRINK 13/14



3906S HEIDBRINK 13/14 SERRATED



3921 HOWARD 2/3

4001S POWER ROOT TIP PICKS SET OF 2



4001 POWER ROOT TIP PICK LEFT DAVIS TEASER 2



4002 POWER ROOT TIP PICK RIGHT DAVIS TEASER 3

4033M POWER FLEX PREMIUM TWIST ROOT PICKS SET OF 3



4033M-1 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR LEFT S/E



4033M-2 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR RIGHT S/E



4033M-3 POWER FLEX TWIST PERIOTOME/ROOT PICK POSTERIOR D/E