

PERIODONTAL SURGERY



C



SURGICAL & KNIVES


Periodontal dental surgery knives are used for a variety of procedures, including creating bevel incisions, contouring the gingiva, excising interproximal tissue, removing excess gum tissue, trimming gum shape, and promoting periodontal tissue regeneration and repair.



1801 K15/16




1802 KGF8



1806 BUCK5/6



1803 KGF11




1804 KSAN1/2



1807 KGF7



1808 KGF9



1805 KUSC19/20



1809 BUCK3/4




1810 KPO1/2



1811 SOLT1/2



1812 KOR1/2



1812L KOR1/2 LONG BLADE



1812M KOR1/2 MICRO BLADE




1813 KCKE



1816D KM3/4

PERIODONTAL FILES


Periodontal files are dental instruments used to break up calculus deposits and other dental work including debridement in deep pockets and scaling and root planing. They can also be used to remove amalgam overhangs and overhanging restorations.



1818 F9/10S



1819 F11/12S



1820 F1/2S



1821 F3/4S

SURGICAL CURETTES

Surgical curettes are used to remove subgingival calculus, plaque, and diseased tissue from the teeth and periodontal pockets during periodontal surgery procedures. They can also be used to remove unwanted debris or cysts from the tooth socket, transfer bone grafts, and clean out infected cavity areas.



1901 PR1/2



1902 KRA1



1903 KRA2




1904 KRA3



1905 KN3/4



1907 KIR4/5




1908 KIR8/9



1918 SS3/6




1919 SS9/10



1920 SS11/12

PERIOTOMES

A perioste is a dental instrument used to remove a tooth by cutting the periodontal ligament (PDL) that holds it in place. The serrated blade cuts the PDL easily, resulting in minimal tissue damage.




1930FSG

X PERIOTOME POWER FLEX & ROOT EXTRACTION SET OF 8



1930-1G POWER FLEX PERIOTOME STR S/E



1930-2G POWER FLEX PERIOTOME CVD S/E



1930-3G X ROOT EXTRACTION TITANIUM SCREW ANTERIOR



1930-4G POWER FLEX PERIOTOME FPT-2 D/E



1930-5G POWER FLEX PERIOTOME FPT-1 D/E



1930-6G POWER FLEX PERIOTOME TWIST D/E




1930-3SG X ROOT EXTRACTION TITANIUM SCREW POSTERIOR



7821G X ROOT EXTRACTION TITANIUM SCREW HANDPIECE ANCHOR 30mm

4032


POWER TWIST COMPLETE SET OF 6



4024 POWER TWIST ELEVATOR STR-R




4023 POWER TWIST ELEVATOR STR-L



3897G POWER TWIST PERIOTOME 4mm POSTERIOR, CVD-R



3896G POWER TWIST PERIOTOME 4mm POSTERIOR, CVD-L



3895G POWER TWIST PERIOTOME 4mm ANTERIOR, STR-R



3894G POWER TWIST PERIOTOME 4mm ANTERIOR, STR-L

SURGICAL CHISELS

Surgical chisels are used to split, remove and reshape bone, reflect flaps, and remove secondary palatal flaps. They can also be used to elevate fractured root tips and remove remaining fibers on teeth as well as slightly remodeling bone margins, removing soft tissue, and for significant deposits of supragingival calculus.



2001 13K/TG 3.5/3mm



2001T 13K/TG 3.5/3mm TITANIUM



2002T 3K/LR 4.5/3.5mm TITANIUM



2002 13K/LR 4.5/3.5mm



2003 KN1/2



2002E 13K/LR 4.5/3.5mm BLACK EDITION



2004 C36/37



2002L 13K/13KL 4/6mm



2005 44/45



2006 TG



2007 SOLT3/4



2008 TG-O



2010 OCH1



2011 OCH2



2014 OCH4



2015 FED11



2016 FED12