

FORCEPS



F


NO FRACTURE



NO TISSUE DAMAGE



EXTRACTION WITH EASE



4974
F1
UPPER
PREMOLARS



4975
F2
LOWER
PREMOLARS



4976
F2
LOWER
PREMOLARS



4977N
F4N
LOWER
ANTERIOR




4979
F6
UPPER
MOLARS



4980C
F7C
LOWER
MOLARS

AMERICAN PATTERN

American pattern dental extraction forceps are typically identified by the shape of connecting point where two forceps beaks connect. They feature hinges that are completely polished out and smooth so that the screw isn't visible.



48011 UPPER CENTRAL



4801AF
1AF
UPPER
CENTRAL



4803
6
LOWER
SPLIT



4804
10H
UPPER
MOLAR



4805
15
LOWER
MOLAR



4806
16S
CHILD
LOWER
MOLAR
COW HORN

EXTRACTION FORCEPS AMERICAN PATTERN



EXTRACTION FORCEPS AMERICAN PATTERN





EXTRACTION FORCEPS ENGLISH PATTERN

ENGLISH PATTERN

English pattern forceps have visible screws in the hinge area. The hinge is vertical, the beaks are at a right angle to the handle and the fulcrum is closer to the working end.



49011 UPPER CENTRAL



EXTRACTION FORCEPS ENGLISH PATTERN

EXTRACTION FORCEPS ENGLISH PATTERN



4955
95
UPPER
MOLAR



4962
168
UPPER
CANINE



4964
233
LOWER
ROOT



4966
MD1
UPPER
WISDOM



4969
MD4
LOWER
MOLAR

EXTRACTION FORCEPS

Extraction forceps are primarily used for the removal of teeth that are decayed, fractured, or impacted. They are also used to remove remaining roots after a tooth has been fractured or previously extracted.



4972 22 L-L/U-R SPECIAL MOLAR



4973
22
L-R/U-L
SPECIAL
MOLAR



4974
F1
UPPER
UNIVERSAL





4974B
PEDO F1
UPPER
UNIVERSAL



4974P
F1P
UPPER
UNIVERSAL



4975
F2
LOWER
UNIVERSAL



EXTRACTION FORCEPS



5103 EXTRACTION TWEEZER ALBRECHT 16cm ROUND END



5104 EXTRACTION TWEEZER ALBRECHT 16cm POINTED END

EXTRACTION FORCEPS

4972S X EXTRACTION FORCEPS SPECIAL MOLAR SET



4972 22 L-L/U-R SPECIAL MOLAR


4973 22 L-R/U-L SPECIAL MOLAR

4994PS X ROOT EXTRACTION FORCEPS SET



4997T 51SR SPECIAL ROOT UPPER

4994T 47SR SPECIAL ROOT LOWER



5033 DR JAMES RIDGE COMPRESSING FORCEP

F-13 FORCEPS

5040
ATRAUMATIC UNIVERSAL EXTRACTION FORCEP SET



4974P
F1P
UPPER
UNIVERSAL



4975P
F2P
UPPER
UNIVERSAL



4976P
F3P
UPPER
ANTERIOR

F-14 FORCEPS



4977P
F4P
LOWER
ANTERIOR

4979P
F6P
UPPER
MOLAR

4980CP
F7CP
LOWER
MOLAR

**5042
FULCRUM FORCEPS KIT**



5042-1 FULCRUM EXTRACTION FORCEP LEFT



5042-2 FULCRUM EXTRACTION FORCEP RIGHT

5042-P1
FULCRUM BLUE
THICK PAD 12PCS



5042-P2
FULCRUM WHITE
THIN PAD 12PCS



4033M-1 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR LEFT S/E



4033M-2 POWER FLEX TWIST PERIOTOME/ROOT PICK ANTERIOR RIGHT S/E



4033M-3 POWER FLEX TWIST PERIOTOME/ROOT PICK POSTERIOR D/E