

IMPLANT & SURGICAL



D

BLACK TITANIUM INSTRUMENTS

Our black titanium instruments are enhanced in strength and durability up to three times more than ordinary tips. They are easy to clean and require less sharpening than other instruments. The titanium coating provides a non-stick surface.



1601T GR1/2



1201BT H6/7 BLACK EDITION HANDLE



1201T H6/7



1203T H5/33



1208BT 204S BLACK EDITION HANDLE



1208T 204S



1209T 204SD




1226T NEVI-2



1322ST MC13/14S POINTED



1415BT COL4R/4L BLACK HANDLE



1415T COL4R/4L



1416BT COL13/14 BLACK HANDLE



1416T COL13/14T



1422T BH5/6



1601BT GR1/2 BLACK HANDLE



1602T GR3/4



1606BT GR11/12 BLACK HANDLE



1606T GR11/12



1607BT GR13/14 BLACK HANDLE



1607T GR13/14



1608T GR15/16



1609T GR17/18



1610T AFG1/2



1614T AFG11/12



1615T AFG13/14



1620T SUB1/2



1621T SUB11/12



1622T SUB13/14

SURGICAL CHISELS

Surgical chisels are used to prepare the bone for dental implants by contouring, smoothing, or removing the bone so that implant material can easily be inserted. They can also be used to split teeth, flatten enamel, and remove secondary palatal flaps.




2001 13K/TG 3.5/3mm



2001T 13K/TG 3.5/3mm TITANIUM



2002L 13K/13KL 4/6mm



2002 13K/LR 4.5/3.5mm



2002T 3K/LR 4.5/3.5mm TITANIUM



2002E 13K/LR 4.5/3.5mm BLACK EDITION



2010 OCH1



2011 OCH2

TITANIUM INSTRUMENTS

Our titanium instruments are enhanced in strength and durability. The titanium coating provides a non-stick surface and an elegant and chic look for any type of scaling, on implants or for calculus.



3506T PERIOSTEALS MOLT9



3509T PERIOSTEALS PR3



3521T PERIOSTEALS HF20



3514T PERIOSTEALS 24G



3526T PERIOSTEALS WOODSON-1 MODIFY



3516T PERIOSTEALS BUSER



3534T PERIOSTEALS PERIO SURGICAL

OSTEOTOMES


An osteotome is a surgical instrument used to cut, compress, or deform bone during dental implant procedures or bone grafting. They can also be used to create implant sites and elevate the sinus floor.



2057 OSTEOTOME 4mm STR (7-10-13-15-18mm)




2058 OSTEOTOME 4mm CVD (7-10-13-15-18mm)



2059 OSTEOTOME 7.5mm STR (7-10-13-15-18mm)



2060 OSTEOTOME 7.5mm CVD (7-10-13-15-18mm)



2061 OSTEOTOME 6.0mm JOVANOVIĆ



2064 RIDGE SPLIT CHISEL 3mm SINGLE/CUT (8-10-13-16mm)



2064D RIDGE SPLIT CHISEL 3mm DOUBLE/CUT (8-10-13-16mm)



2065 RIDGE SPLIT CHISEL 6mm SINGLE/CUT (8-10-13-16mm)




2065D RIDGE SPLIT CHISEL 6mm DOUBLE/CUT (8-10-13-16mm)




2066 RIDGE SPLIT CHISEL 9mm SINGLE/CUT (8-10-13-16mm)



2066D RIDGE SPLIT CHISEL 9mm DOUBLE/CUT (8-10-13-16mm)



2067
FORCEP BONE BLOCK
19cm STR



2067N
FORCEP BONE BLOCK
19cm CVD



2068
FORCEP CORTICAL BONE
BLOCK 19cm CVD

SINUS LIFT INSTRUMENTS

Sinus lift instruments are utilized for gently lifting the sinus membrane away from the floor and lateral portion of the maxillary sinus. Coming in a variety of shapes and sizes, this ensures clean lift at any angle necessary.



2070S SINUS LIFT MEMBRANE ELEVATION KIT SET OF 7 PCS



1970 SINUS BONE GRAFT 3.87SP/2.8mm (2-4-5-6-8-10mm)

2062 SINUS LIFT 6

2063 SINUS LIFT 7

2069 MUSHROOM 2.8mm/3.8mm (2-4-5-6-8-10)

2035T BONE CARRIER TITANIUM

2038T BONE GRAFT/CARRIER 2-4-6-8-10mm TITANIUM

7503S CASSETTE, DOUBLE HINGE 6 HOLD SINGLE TALL RACK 90X203X31mm



COMMONLY USED IMPLANT INSTRUMENTS

Aid your implant procedures with the most commonly used implant instruments. Having a complementary set of instrumentation that focuses on all aspects of implant placement is necessary for success in any placement. Bone graft carriers, pluggers, tunneling periostomes, and more, every instrument needed in implant protocols is included in this section of our catalog.



1923T PT2 TITANIUM



1924T PT1 TITANIUM



2033A BONE GRAFT PACKER
1.0/1.4mm (8-11-13-15-18mm) TITANIUM



1928 MS1



2033B BONE GRAFT PACKER
1.7/2.0mm (8-11-13-15-18mm) TITANIUM



2032 BONE PLUGGER IMPLANT
LABIANCA 3X1.5mm



2033C BONE GRAFT PACKER
2.0/2.7mm (8-11-13-15-18mm) TITANIUM



2033 BONE PLUGGER IMPLANT
JOVANOVIC 3mm



2033N BONE GRAFT PACKER
3.3/4mm (8-11-13-15-18mm)




3533N PERIOSTEAL P25XL 6EZ



2033NT BONE GRAFT PACKER
3.3/4mm TITANIUM



2036S BONE GRAFT PACKER
4/6mm SERRATED



2033NTS BONE GRAFT PACKER
3.3/4mm TITANIUM SERRATED



2038T BONE GRAFT/CARRIER
2-4-6-8-10mm TITANIUM



2034S BONE GRAFT PACKER
3/5mm SERRATED



2043T PPAEL TITANIUM



2035T BONE CARRIER TITANIUM



2044T PPAELA TITANIUM



1925T
PT3
TITANIUM



1926T
PT4
TITANIUM



1927T
PT5
TITANIUM



2037T
BONE GAUGE
4-6-8-10-12-14-16mm
TITANIUM

2030 BONE SPOON IMPLANT BSJO5




2071 KERRISON FORCEP 3mm
2072 KERRISON FORCEP 4mm




2087 BONE CRUSHER

2083
BONE CRUSHER
(MORSELIZER)



1941 BONE SYRINGE 5mm
2079 BONE SYRINGE 3mm
2080 BONE SYRINGE 3.5mm
2081 BONE SYRINGE 4mm
2082 BONE SYRINGE 6mm



2073 MICRO SURGERY FORCEP 75mm
2074 MICRO SURGERY FORCEP 95mm
2075 MICRO SURGERY FORCEP 110mm



2091 BONE COLLECTOR



2093 BONE COLLECTOR CHISEL STR



2093B BONE COLLECTOR STR BLADE ONLY



2093A BONE COLLECTOR CHISEL CVD



2093AB BONE COLLECTOR CVD BLADE ONLY



2092 BONE COLLECTOR WITH FILTER



2092F BONE COLLECTOR FILTER ONLY



2094
BONE MIXING BOWL





4613 CASTRO CALIPER STR 8.5cm



4613C CASTRO CALIPER CVD 10cm



4613L CASTRO CALIPER CLC 17cm STR



4613LA CASTRO CALIPER CLC 17cm CVD



4618 RULER 15cm 6in



4613RA CASTRO CALIPER WEISS MODIFIED CVD CLC 17cm CONTACT FRONT



4614 BONE RIDGE MAPPING CALIPER 0-35mm



4625 BONE RIDGE MAPPING CALIPER IMPLANT



4628 MICRO BOLEY CALIPER GAUGE



4411 BONE RONGEURS MICRO FRIEDMAN 15 14.5cm



4412 BONE RONGEURS MICRO FRIEDMAN 90 14.5cm



4413 BONE RONGEURS MICRO FRIEDMAN 90 14.5cm



6436 WEIDER TONGUE RETRACTOR 30mm



6436L WEIDER TONGUE RETRACTOR 35mm



6445 LIP RETRACTOR

7425 BONE TRAY SQUARE 65x58x13mm



6442 CHALAZION 9mm



6443 CHALAZION 24x14mm



6443S CHALAZION 16x10.5mm



6444 CHALAZION 24x14mm



7426 BONE TRAY BOWL 58mm DIAx28mm DEPTH



5309 PEAN
BABY CRILLE
STR 14cm
H-1



5309C PEAN
BABY CRILLE
CVD 14cm
H-2



5926 CASTRO
STR 10cm




5926C CASTRO
CVD 10cm



5929 CASTRO STR 14cm



5929C CASTRO CVD 14cm



5944 CASTRO MICRO 45° ANGLE 18cm

TISSUE PLIERS

Manage your tissue effectively and precisely with our Tissue Pliers. Every type of tissue plier, for all types of soft tissue management.



6501 ADSON STR 12cm



6501A ADSON CVD TEETH 12cm



6501L ADSON STR TEETH 15cm



6502 ADSON STR TEETH 1x2 12cm



6502A ADSON CVD TEETH 1x2 12cm



6502AL ADSON CVD TEETH 1x2 15cm



6502B ADSON BROWN 12cm



6502BL ADSON BROWN 15cm



6503 SEMKEN-TAYLORS STR REGULAR 12.5cm



6505 SEMKEN-TAYLORS STR TEETH 1X2 12.5cm



6505L MICROSURGICAL STR TEETH 1X2 15cm



6506 SEMKEN-TAYLORS CVD TEETH 1X2 12.5cm



6506L MICROSURGICAL CVD TEETH 1X2 15cm



6507 FORCEPS POINT 12cm



6508 FORCEPS POINT 11.5cm



6509 FORCEPS POINT CVD 11.5cm



6510 FORCEPS POINT CVD 11.5cm



6511 SUTURE CORN TISSUE FORCEPS CVD 16cm



6511E SUTURE CORN TISSUE FORCEPS ATRAUMATIC CVD 18cm
BLACK EDITION



6511N SUTURE CORN TISSUE FORCEPS MICRO SUTURE CVD 16cm



6511S SUTURE CORN TISSUE FORCEPS STR 16cm



6512 GOLDMAN-FOX 15cm L



6513 GOLDMAN-FOX 15cm R



6515 SELF-LOCKING 15cm



6517 LONDON-COLLEGE



6517G LONDON-COLLEGE TITANIUM



6517L LONDON-COLLEGE 15cm



6518 MERIAM



6519 PERRY 13cm



6520 MERIAM SELF-LOCKING 16cm



6523 SELF-LOCKING




6525 GERALD STR 18cm



6525TC GERALD CVD 18cm MICRO



6526 TOMOR 16cm



6526L TOMOR 20cm



6527 TP16



6528 TP50L



6531 DEBAKEY ATRAUMATIC STR 18cm



6531B DEBAKEY ATRAUMATIC STR 18cm BLACK TITANIUM



6531C DEBAKEY ATRAUMATIC CVD 18cm



6531CB DEBAKEY ATRAUMATIC CVD 18cm BLACK TITANIUM



6532 GERALD STR REGULAR 18cm



6532C GERALD CVD REGULAR 17.5cm



6533 SPECIAL MEMBRANE GRAFTING



6534 ABUTMENT HOLDER
HOLE SIZE 1-2.5mm/1.6-3.2mm



SECTIONAL FORCEPS



4629
FINGER RULER 38mm



6529
ALLISON TFB
12.5cm



6530
ALLISON TFS
15cm

7412
ENDODONTIC BOX
10cmX4.5cmX4.2cm



5623
NERVE CANAL
NARROW
12.5cm



5106E ROOT FRAGMENTS 45°



5107E ROOT FRAGMENTS 90°

TISSUE PUNCHES

Remove soft tissue with precision and improve accuracy and access in implant placement with our Dental Tissue Punches.



2076 STR 4mm



2077 STR 5mm



2078 STR 6mm



1955A CVD 2mm



1956A CVD 3mm



2076A CVD 4mm



2077A CVD 5mm



2078A CVD 6mm

ACCU BERNARD SPADE SET

Our Accu Bernard Spade Set is manufactured with surgical grade stainless steel coated with blue titanium. The fine micro serrations easily loosen the tooth from the socket and periodontal ligaments. The syndesmotomes can be used to elevate small roots.



3771 ACCU LUX-TOME SPADE




3790 ACCU BERNARD SPADE OFFSET IN



3791 ACCU BERNARD SPADE OFFSET OUT



3792 ACCU BERNARD SPADE LEFT




3793 ACCU BERNARD SPADE RIGHT

PERIOSTEALS

Periosteals are used to separate and elevate the periosteum after removing gingival tissue. They can also be used for tooth extraction, bone grafting and sinus lifts.

D-25 IMPLANT & SURGICAL



D-26 IMPLANT & SURGICAL



3520 HOURIGAN2



3522 MOLT5



3523 SELDIN23



3523N SELDIN23 BLUNT



3525 OHL10



3531 P25



3531N P25 NEW HANDLE



3532 P25L



3533 P25XL



3533E P25XL BLACK EDITION



3533N P25XL 6EZ



3540 PRICHARD DOUBLE END TITANIUM

SURGICAL CURETTES (SPOON)


Surgical curettes are used for curettage and cyst removal as well as craping or debriding tissue. The spoon tip allows for precision in removing unwanted or unhealthy tissue without damaging surrounding areas.



4101 LUCAS84



4101S LUCAS84S SERRATED




4101SE LUCAS84S SERRATED BLACK EDITION



4102 LUCAS85



4102S LUCAS85 SERRATED




4102SE LUCAS85S SERRATED BLACK EDITION



4103 LUCAS86




4103AS LUCAS86AS SERRATED OFF ANGLE



4103ASE LUCAS86AS SERRATED OFF ANGLE BLACK EDITION



4103S LUCAS86S SERRATED




4103SE LUCAS86S SERRATED BLACK EDITION



4104 MILLER8



4105 MILLER9




4106 MILLER10



4107 LUCAS87



4108 LUCAS88



4109 MILLER11



4110 MILLER12



4111 MOLT2/4



4112 MOLT2



4113 MOLT4



4114 MOLT5L



4115 MOLT5R



7501 CASSETTE 5 HOLD

4126
POWER SERRATED CURETTE KIT




4101S LUCAS84S SERRATED



4102S LUCAS85 SERRATED



4103AS LUCAS86AS SERRATED
OFF ANGLE



4103S LUCAS86S SERRATED

4126E
POWER SERRATED CURETTE KIT BLACK EDITION




4101SE LUCAS84S SERRATED BLACK EDITION



4102SE LUCAS85S SERRATED BLACK EDITION



4103SE LUCAS86S SERRATED BLACK EDITION



4103ASE LUCAS86AS SERRATED OFF ANGLE BLACK EDITION


BONE CHISELS

Bone chisels are used to shape and smooth bone tissue around teeth. They are also used to remove excess bone, trim alveolar bone, split bone, and prepare the bone for dental implants.



BONE FILES

Bone files are used to remove small bone segments, tooth decay, and shape dental fillings. They can also be used to smooth bone and narrow the canals of the interior of a tooth. Their serrated edges make it easy to maneuver between soft tissue and bone.



MALLETS

Mallets are used to impact dental materials and instruments during dental procedures. They are also used to split teeth, remove bone, and break up bone tissue. They are often used with chisels or wedges.



4223N MEAD1



4223M MEAD MINI



4226 MEAD MINI
FOR IMPLANT



4223MT MEAD MINI REPLACEMENT TIPS (2)
4223T MEAD REPLACEMENT TIPS (2)

CHISELS & GOUGES

Chisels are used to remove, smooth, and contour bone so that the implant can be placed. Gouges are used to loosen teeth from the gums to make extraction easier and to prevent damage to surrounding tissue and teeth.



4503 FREER 16cm 4mm CVD



4504 LUCAS 16cm 4mm



4507 14cm 6mm



4508 PARTCSH 14cm 5mm



4510 PARTCSH 14cm 4mm

BONE RONGEURS

Bone rongeurs are used to cut or remove small pieces of bone or tissue. In addition to being used to trim and contour bone they can also be used to open a window in the bone to access tissue underneath.



4401
BLUMENTHAL
15cm



4402
BLUMENTHAL
45° 14.5cm



4403
BLUMENTHAL
90° 14.5cm



4405
BEYER
DOUBLE
ACTION
18cm



4406
BEYER
DOUBLE
ACTION
RUSKIN-MINI
15.5cm



4408
FRIEDMAN
14cm



4409
MINI
FRIEDMAN 15
12cm



4411
MICRO
FRIEDMAN 15
14.5cm



4411L
MICRO
FRIEDMAN 15
15.5cm



4412
MICRO
FRIEDMAN 90
14.5cm



4413
MICRO
FRIEDMAN 45
14.5cm

ROOT SPLINTER

Root splinters are used to extract teeth and roots, and can be used on both the upper and lower gums. They are designed to prevent fractures and are effective on roots due to their serrated tips, curved jaws, and slim profile.



5105 ROOT FRAGMENTS STR



5106
ROOT
FRAGMENTS
CVD 45°



5106E
ROOT
FRAGMENTS
45°



5107
ROOT
FRAGMENTS
CVD 90°



5107E
ROOT
FRAGMENTS
90°



5108
PEET
SHARP
CVD 15°



5109
PEET
SHARP
CVD 45°



5110
PEET
SHARP
CVD 90°

HEMOSTATIC FORCEPS

Hemostatic forceps are used to clamp blood vessels, control bleeding, and move gauze and other materials around the mouth. They can also be used to hold and move small objects like bone and tooth fragments.



5301C
BABY
MOSQUITO
CVD 10cm

5301
BABY MOSQUITO STR 10cm



5302C
MICRO-
HALSTED
CVD 12.5cm

5302
MICRO-HALSTED STR 12.5cm



5303C
HALSTED-
MOSQUITO
CVD 12.5cm

5303
HALSTED-MOSQUITO STR 12.5cm



5305C
MICRO-
HALSTED
CVD 12.5cm

5305
HALSTED-MOSQUITO STR 12.5cm



5307
PEAN 12.5cm



5309C
PEAN
BABY-CRILE
CVD 14cm H-2

5309
PEAN BABY-CRILE STR 14cm H-1



5310C
CRILE/
KELLY
CVD 14cm

5310
CRILE/KELLY STR 14cm



5311C
BABY CRILE/
KELLY KOCHER
CVD 14cm

5311
BABY CRILE/KELLY KOCHER
STR 14cm



5312C
CHILE
KOCHER
CVD 16cm

5312
CHILE KOCHER STR 16cm



5313C
HALSTED
CVD 18.5cm

5313
HALSTED STR 18.5cm



5315C
HALSTED
CVD 18.5cm

5315
HALSTED STR 18.5cm



5316
BABY-MIXTER CVD 18cm



5318
IMPLANT ABUTMENT
HOLDER CVD 14.5cm

MATHIEU PLIERS

Mathieu pliers are used to hold, place, and tie elastics and ligatures, and to guide needles during suturing. are lightweight and flexible, with a spring function that allows the user to lock and unlock the jaws.



5407 MATHIEU
14cm



5407F MATHIEU
FINE 14cm



5407U MATHIEU
US ORTHO 14cm



5414 MATHIEU
DELICATE STR 15cm



5415 MATHIEU
DELICATE CVD 15cm



5418 MATHIEU
14cm



5419 MATHIEU
THIN 14cm



5501N MATHIEU T/C THIN
BLACK TITANIUM 16cm



5502 MATHIEU
T/C WIDE 14cm



5502L MATHIEU
T/C WIDE 17cm

NEEDLE HOLDERS

Needle holders are used to grasp and manipulate needles during suturing. Needle holders can also be used in other dental procedures where precise needle control is required.



5402 CRILE-MURRAY 15cm



5403
MAYO-
HEGAR
15cm



5403L
MAYO-
HEGAR
16cm



5403W
MAYO-
HEGAR
WIDE
15cm



5408 BOYNTON 12cm



5409C CASTROVIEJO CVD T/C 16cm



5409B CASTROVIEJO STR
BLACK TITANIUM 16cm



5416 CASTROVIEJO
MATRIX CVD 13cm



5505L
CRILE-WOOD T/C 18cm



5505TN
CRILE-WOOD STR T/C 15cm



5506
CRILE-WOOD CVD T/C 15cm



5508
RYDER STR T/C 18cm



5509
DERF STR T/C 11.5cm



5512
HALSEY STR T/C 13cm



5521
MAYO-HEGAR STR T/C 18cm



5521L
MAYO-HEGAR STR T/C 20cm



5545
CRILE-WOOD "S" CVD HANDLES T/C 15cm



5504B OLSEN-HEGAR
BLACK TITANIUM STR 15cm



5505B CRILE-WOOD
BLACK EDITION STR 15cm



5505BN CRILE-WOOD
BLACK TITANIUM STR 15cm



5505LB CRILE-WOOD
BLACK EDITION STR T/C 18cm



5505TB CRILE-WOOD
BLACK TITANIUM
STR T/C 15cm



5505TNB CRILE-WOOD
BLACK EDITION
STR T/C 15cm



5506B CRILE-WOOD
BLACK EDITION
CVD T/C 15cm



5509B DERF
BLACK EDITION
STR T/C 11.5cm



5512B HALSEY
BLACK EDITION
STR T/C 13cm



5520B MAYO-HEGAR
BLACK EDITION
STR T/C 16cm



5521B MAYO-HEGAR
BLACK TITANIUM
STR T/C 18cm



5513B
CASTROVIEJO
BLACK EDITION
STR 16cm



5513BN
CASTROVIEJO
BLACK TITANIUM
STR 16cm



5514B
CASTROVIEJO
BLACK EDITION
STR 18cm



5514BN
CASTROVIEJO
BLACK TITANIUM
STR 18cm



5515B
CASTROVIEJO
BLACK EDITION
CVD 16cm



5515BN
CASTROVIEJO
BLACK TITANIUM
CVD 16cm



5516B
CASTROVIEJO
BLACK EDITION
CVD 18cm



5516BN
CASTROVIEJO
BLACK TITANIUM
CVD 18cm



5527
DOUBLE
ACTION JAWS
14cm



5531
CASTROVIEJO
SPECIAL
CVD 45°
16.5cm



5531N
MATRIX
PLACEMENT &
REMOVAL
CVD

GUM SCISSORS

Gum scissors are used primarily for cutting soft tissue in the oral cavity, such as gingival tissue. They are ommonly used in periodontal surgeries, such as gingivectomy or flap surgeries, to trim excess gum tissue.



5904 KELLY STR 16cm

5905
KELLY
CVD
16cm



5906 QUINBY STR 12.5cm

5906C
QUINBY
CVD
12.5cm



5907 DEAN 17.5cm



5907N DEAN
GOLD TITANIUM
SCREW 17.5cm



5908C
METZENBAUM
CVD
14.5cm

5908 METZENBAUM
STR 14.5cm



5909 LOCKLIN 16cm



5911C
GOLDMAN-FOX
CVD
12.5cm

5911 GOLDMAN-FOX
STR 12.5cm



5916 GOLDMAN-FOX
SPECIAL CVD 12.5cm



5917 LA GRANGE
CVD 11.5cm



5919 IRIS STR 11.5cm



5920
IRIS
CVD
11.5cm



5920F IRIS MICRO
CVD 11.5cm



5922 SUTURE 15cm



5905B
KELLY
SUPER CUT
BLACK EDITION
CVD 16cm

5904B KELLY SUPER CUT
BLACK EDITION STR 16cm



5907B DEAN SUPER CUT
BLACK EDITION 17.5cm




5907BN DEAN
POWER BOW 17.5cm



5908B METZENBAUM
BLACK EDITION STR 14.5cm



5908CB
METZENBAUM
BLACK EDITION
CVD 14cm



5911B GOLDMAN-FOX
BLACK EDITION STR 12.5cm

5911CB
GOLDMAN-FOX
BLACK EDITION
CVD 12.5cm



5911BN GOLDMAN-FOX
POWER BOW STR 12.5cm

5911CBN
GOLDMAN-FOX
POWER BOW
CVD 12.5cm



5916B GOLDMAN-FOX
SPECIAL CVD
BLACK EDITION 12.5cm



5917B LA GRANGE
SUPER CUT
BLACK EDITION CVD 11.5cm




5917BN LA GRANGE
POWER BOW CVD 11.5cm



5919BN IRIS
POWER BOW STR 11.5cm




5920BN
IRIS
POWER BOW
CVD 11.5cm



5919B IRIS SUPER CUT
BLACK EDITION STR 11.5cm



5920B
IRIS
SUPER CUT
BLACK EDITION
CVD 11.5cm



5922B SUTURE SUPER CUT
BLACK EDITION 15cm



5904S KELLY SUPER CUT
STR 16cm



5907S DEAN
SUPER CUT 17.5cm



5908S METZENBAUM
SUPER CUT STR 14.5cm



5908SP METZENBAUM
POINT FRONT STR 14.5cm



5911S GOLDMAN-FOX
SUPER CUT STR 12.5cm



5916S GOLDMAN-FOX
SUPER CUT TUNGSTEN
BLACK EDITION 12.5cm



5919S IRIS
SUPER CUT STR 11.5cm



5920S IRIS
SUPER CUT
CVD 11.5cm

5908SC
METZENBAUM
SUPER CUT
CVD 14.5cm

5908SCP
METZENBAUM
POINT FRONT
CVD 14.5cm

5911S
GOLDMAN-FOX
SUPER CUT
CVD 12.5cm



5923
SUTURE SLIM
12cm



5924
CROWN
STR 12cm



5926
CASTRO
STR 10cm



5926B
CASTRO
BLACK
TITANIUM
STR 10cm



5926C
CASTRO
CVD 10cm



5929
CASTRO
STR 14cm



5929C
CASTRO
CVD 14cm



5929CS
CASTRO
MICRO
CVD 14cm



5929S
CASTRO
MICRO
STR 14cm




5940
CASTRO
MICRO
LOCKLIN
45° 18cm



5944
CASTRO
MICRO
OFF ANGLE
45° 18cm

TUNGSTEN CARBIDE INSERT SCISSORS


Harder than traditional stainless steel, cut through the toughest of materials with our Tungsten Carbide Scissors.



6001 METZENBAUM
ROUND T/C STR 14.5cm

6002

METZENBAUM
ROUND
T/C CVD 14.5cm



6003 METZENBAUM
POINT T/C STR 14.5cm

6004

METZENBAUM
POINT
T/C CVD 14.5cm



6003BN METZENBAUM
POWER BOW
T/C STR 15cm

6004BN

METZENBAUM
POWER BOW
T/C CVD 14.5cm



6008 WIRE CUTTING
T/C 12cm



6009 IRIS STR
T/C 11.5cm

6010

IRIS CVD
T/C 11.5cm



6011 DEAN
T/C 17.5cm



6012 LOCKLIN
T/C 16.5cm



6014 KELLEY STR
T/C 16cm

6015

KELLEY CVD
T/C 16cm

SKIN HOOKS

Skin hooks are used hold back the edges of skin during dental procedures. They can also be used to elevate tissue and mobilize it as well as provide gentle horizontal traction.



6312 GILLIES 18cm



6345 MICRO SURGERY 28K-L



6346 MICRO SURGERY 28K-R

RETRACTORS

Retractors are used to hold back soft tissues, such as gums and cheeks, to provide better visibility and access to the treatment area. Mouth gags are used to prop the mouth open during procedures.

6342 FARABEUF 12.5cm

6343 FARABEUF 15cm

6344 ROUX 16.5cm

6420 UNIVERSITY OF MINNESOTA

6420B UNIVERSITY OF MINNESOTA
CONCAVE FRONT

6420R UNIVERSITY OF MINNESOTA
ROUND FRONT

6420A UNIVERSITY OF MINNESOTA
FLAT WITH SERRATION

6420N UNIVERSITY OF MINNESOTA
SPECIAL HANDLE HEAVY HOLD

6421 ALM 7cm

6423 CAWOOD MINNESOTA

6428 BISHOP

6426 MOUTH GAG ADULT
6426T BLUE REPLACEMENT TIPS (ADULT SIZE)
6427 MOUTH GAG CHILD
6427S MOUTH GAG PEDO
6427T BLUE REPLACEMENT TIPS (CHILD SIZE)

6424 COLUMBIA 14cm

6425 ORRINGER HANDS FREE

6429 HUMAN ADULT

6430 HUMAN CHILD

6433 VESTIBULUM LIP LARGE

6434 VESTIBULUM LIP SMALL

6434N UNIVERSAL CHEEK
AND LIP SMALL

6435 BAYONET 15cm



6438 SPANDEX
SMALL 10cm



6440 S/E 9cm



6446
MOUTH PROPS
ADULT



6447
MOUTH PROPS
CHILD



6445 LIP RETRACTOR



6448 AUSTIN



6449 STERNBERG RETRACTOR



6448A AUSTIN DOUBLE



6450 CHEEK RETRACTOR



6451 CHEEK DOUBLE SIDE



6452 SINUS SOFT TISSUE HOLD AND EXPANDER 10mm



6453 SINUS SOFT TISSUE HOLD AND EXPANDER 15mm



6454 CHEEK AND TONGUE SPOON TYPE

SCALPEL HANDLES

Scalpel handles are used in conjunction with mess blades for making precise incisions in soft tissue. Our handles are designed for comfort and control during procedures.



6901 4 FLAT



6902 3G FLAT



6904 5A ROUND CVD



6905 5A ROUND STR



6905D 5A ROUND DOUBLE HEAD



6906 7 THIN HANDLE



6907 3L FLAT LONG



6920 6EZ



6920A PALATAL HARVESTING 1.5mm GAP/1mm



6920B PALATAL HARVESTING 1.5mm GAP/2mm



6920BE IMPLANT SURGERY BLACK EDITION



6920G 6EZ DR. GEORGE TIP SHAPE 45° CUT



6924 S/E 360° BLADE ANGLE TURN AROUND



6934 BLACK EDITION FOR MICRO BLADE



6935 D/E 45° ORBAN TYPE OFF ANGLE
S/E 360° BLADE ANGLE TURN AROUND

6920S
PALATAL HARVESTING KIT




6920A SCALPEL HANDLE



6920B SCALPEL HANDLE

6966 EASY BLADE REMOVAL FORCEP



6923
SCALPEL HANDLE
180° ADJUSTABLE
16.2cm



SYRINGES

Syringes are used for injecting anesthetics from a cartridge. Three-way syringes supply air, water or a mist through separate interval channels created by combining pressurized air with water flow.



7001 INTRALIGMENTARY 1.8ml



7003 1.8ml WITH HARPOON



7017 PAROJECT



7018L ASPIRATED 2.2ml WING




7018LN ASPIRATED 2.2ml WING



7018N ASPIRATED 1.8ml WING



7019 3-WAY



7019T 3-WAY TIPS (2)